



CP9xx

Mustervisualisierung 15" und 24"



CP9xx – Mustervisualisierung


Stil: Dark




Die Bedienoberfläche der COMTRAXX® CP9xx Tableaus kann individuell an die jeweiligen Anforderungen und das Corporate Design angepasst werden.

Hinweis:
Die unterschiedlichen Stile sind frei wählbar.



Stil: 3D-Dark



Stil: White



CP9xx – Mustervisualisierung





CP9xx – Mustervisualisierung



The screenshot displays the CP9xx control interface. On the left is a vertical navigation menu with icons for Home, Lighting, AHU (highlighted), Timer, RGB Lights, and IPS. The main area is divided into two sections: 'Ventilation' and 'Air Conditioning'. Under 'Ventilation', there are two toggle switches: 'Full UCV Mode' (ON) and 'Conventional Mode' (OFF). The 'Air Conditioning' section shows a temperature display with '0°C' and a vertical slider. Below this, three data boxes are visible: 'Pressure OR' (0bar), 'Temperature Corridor' (21,4°C), and 'Pressure Corridor' (0bar). A green status bar at the bottom indicates 'System OK'. At the very bottom, a row of system status indicators shows 'Ok' for 'IPS System', 'UPS System', 'AC / Vent', and 'Med. Gases'.

This screenshot shows the CP9xx control interface for the 'Timer' function. The left navigation menu is identical to the previous screenshot, with 'Timer' highlighted. The main area features two large digital displays: 'Timer 1' showing '00:0230' in green and 'Elapsed Timer' showing '02:1500' in orange. To the right of each display are three buttons: 'Start', 'Stop', and 'Reset'. At the bottom of the timer section, there are three control buttons: '+/- 1h', '+/- 1min', and '+/- 1sec'. A yellow banner in the top right corner reads 'Optional: Auf Anfrage'. The green 'System OK' status bar and the bottom system status indicators are also present.

CP9xx – Mustervisualisierung



CP9xx – Mustervisualisierung



CP9xx – Mustervisualisierung

Lamp 1

Control all	Brightness	Light area	Color
On	+	Wide	3700 K
Off	85%	Narrow	4400 K
Sync	-	Ambilite	5100 K

Lamp 2

Control all	Brightness	Light area	Color
On	+	Wide	3700 K
Off	65%	Narrow	4400 K
Sync	-	Ambilite	5100 K

Lamp 3

Control all	Brightness	Light area	Color
On	+	Wide	3700 K
Off	85%	Narrow	4400 K
Sync	-	Ambilite	5100 K

Lamp

Camera

System OK

Camera 1

On	Pause	Zoom	Focus	Iris
Off		+	+	+
	Rotate	-	-	-

Camera 2

On	Pause	Zoom	Focus	Iris
Off		+	+	+
	Rotate	-	-	-

Lamp


Camera

System OK



Bender GmbH & Co. KG

Londorfer Straße 65 • 35305 Grünberg • Germany
Tel.: +49 6401 807-0 • info@bender.de • www.bender.de



BENDER Group